

TRANSPARENCY

IN CUSTOM PROCEDURE TRAY:
A STANDARD REQUIREMENT

“Make transparency become a standard requirement” is the principle that drives every business choice in STS. A real ethical and professional approach that has been expressed in adopting, for the production and certification of the Procedural Sets, article 12 of the directive n. 93/42/EEC; instead of article 11 commonly adopted by the most of the manufacturers operating in the European market.



**WE ANTICIPATE
THE FUTURE**

WATCH OUR VIDEOS ON OUR
YOUTUBE CHANNEL



STS Medical Group

Société à responsabilité limitée
4, rue du Fort Wallis / L-2714 Luxembourg
Phone +352 27 36 55 87 / Fax : +352 281152 299
info@stsmedicalgroup.com / www.stsmedicalgroup.com

The Bulgarian Manufacturing Hub

Industrial Area Sokolovetz - 2800 Sandanski - Bulgaria
Phone: +359 (0) 746 31555 - Fax: +359 (0) 746 38798
www.sts-bmh.com

Ditta Luigi Salvadori S.p.A.

Viuzzo del Piscetto 6-8 - 50018 Scandicci - Italy
Tel: +39 (0) 55 75441 - Fax: +39 (0) 55 753686
www.luigisalvadori.it

Sengewald Klinikprodukte GmbH

Adlerstraße 2 - 83101 Rohrdorf - Germany
Tel: +49 (0) 8031 7224 0 - Fax: +49 (0) 8031 7224 295
www.sengewald.de

MSP Schmeiser GmbH

Siemensstraße 14 - 72160 Horb - Germany
Tel: +49 (0) 7451 5535 0



STS MEDICAL GROUP
SECUTRAY®
EFFICIENCY & CUSTOMIZATION
TRANSPARENCY
IN CUSTOM PROCEDURE TRAY
A STANDARD REQUIREMENT



- quantity**
Number of items in each pack
- article code**
Unique STS material identification number
- description**
Detailed description of the component.
- manufacturer**
Original manufacturer of the component. Gives you the guarantee of the quality.
- CE mark**
Allows to identify the notified body if involved.
- risk class**
Identify the level of risk which the component belongs.
- country of origin**
Highlight the actual country where the product is manufactured.

**TRANSPARENCY STARTS FROM
 STS CPTS LABEL**

Secu-Tray is STS's product line for Custom Procedure Trays (CPTs). In advance of EU Regulation, STS has decided to upgrade the quality and quantity of information listed on the labels of its products to provide the highest level of transparency in the CPT market.

In fact, Secu-Tray CPT labels offer the highest level of transparency for each component of any CPT in the European market.

On Secu-Tray labels, the user can immediately see product description, article code, original manufacturer, CE mark, notified body (if required), risk class and country of origin.

Providing detailed and high level information increase transparency and allow the customer to be in control, in fact replacing items without explicit customer approval is not possible.

The Secu-Tray CPT label ensures that the customer is always in control of all medical devices in a CPT, thereby creating a higher-quality surgical procedure and enhancing patient protection.

STS MEDICAL GROUP®
SECUTRAY
 EFFICIENCY & CUSTOMIZATION
Ivom Op Set
1X REF 662 29295
 Treatment Unit according to Article 12 of the Medical Device Directive 93/42/EEC

Q.TY	ARTICLE CODE	DESCRIPTION	MANUFACTURER	CE MARK	RISK CLASS	COUNTRY OF ORIGIN
1	647580-9	Drape 75 x 80 cm, Mediwell 2-ply unfolded	Sengewald	CE	I	BG
2	175015	Copolymer examination glove, Size L	Meditrade	0123	Is	CN
1	575403	Tray, 12 x 14 x 2,5 cm, blue	DeRoyal	CE	I	US
1	9007610	Measuring and fixation instrument 3,5/4 mm, blue	Scilo	CE	I	CN
2	325500	Gauze compresses 10 x 10 cm, 8-ply	Meditrade	CE	I	CN
2	900560	Ocular sticks, Ø 5 mm, 66 mm	Lohman&Rauscher	0123	Ila	DE
2	750300	Gauze balls, plum-size	Meditrade	CE	I	CN
1	900129	Universal eye shield, transparent	Scilo	CE	I	CN
1	9002100-2	Eye speculum, open valves	Miqraz	CE	I	PK
1	555771-9	Eye drape 40 x 50 cm with perforated incision film	Sengewald	CE	I	BG

Ivom Op Set
 REF 662 29295 LOT 123456 2022-04-30
 (01)04034393285283(10)123456(17)220430

Ivom Op Set
 REF 662 29295 LOT 66216138-01 2022-04-30

Ivom Op Set
 REF 662 29295 LOT 66216138-01 2022-04-30

Ivom Op Set
 REF 662 29295 LOT 123456 2022-04-30

Ivom Op Set
 REF 662 29295 LOT 123456 2022-04-30

LOT 123456 YYY-MM 2022-04-30 STERILE EO LATEX

STS MEDICAL GROUP®
 IMPROVING HEALTHCARE THROUGH EFFICIENCY

PROCEDURE PACK
 VALIDATION
 ASSURANCE
 STERILIZATION

Sengewald Klinikprodukte GmbH
 Adlerstrasse 2
 83101 Rohrdorf-Thansau (Germany)

www.stsmedicalgroup.com
 MADE IN GERMANY

TRANSPARENCY BENEFITS

- Components** are clearly CE marked and the **original manufacturer is visible on the label.**
- No silent changes of components** from set packer side are possible.
- Regulatory safety in the light of new EU Medical Device Regulations (MDR) which came into effect last May 2021 as **only certified medical devices are used in sets.**
- What you choose is what you get:** you are in control!